



110 West Cliff Avenue  
Spokane, WA 99204

509.755.8600  
800.541.7891

# TEST CHANGE ALERT #363

---

**July 01, 2010**

## Summary Of Changes

TestCode(s)	Test Description
GLBMAB .....	GLOMERULAR BASEMENT MEMBRANE AB (Reference Range & Notes)



110 West Cliff Avenue  
Spokane, WA 99204

509.755.8600  
800.541.7891

# TEST CHANGE ALERT #363

## July 01, 2010

The following tables reflect revisions only; other existing data remain unchanged.

### GLBMAB

### GLBMAB

### GLOMERULAR BASEMENT MEMBRANE AB

(Reference Range & Notes)

order code

flexilab code

Effective	07/02/10			
Reference Ranges	<i>Glomerular Basement Membrane Ab</i>	<i>Negative Weak Positive Mod-Strong Pos</i>	<i>0-20 21-30 GT 30</i>	<i>Units</i>
	<p><i>This test is designed for the in-vitro measurement of specific IgG autoantibodies against the glomerular basement membrane (GBM). It is intended as an aid in the diagnosis of Goodpasture's syndrome. Some patients with other renal diseases may exhibit positive results. Glomerular Basement Membrane antibodies are not found in normal healthy individuals.</i></p> <p><i>Results were obtained with the QUANTA Lite GBM ELISA assay. Values obtained from different manufacturers' assays cannot be used interchangeably. The magnitude of the reported IgG levels cannot be correlated to an endpoint titer.</i></p>			

[PAML Web Test Directory](#)